

---

# Infant Circumcision Training Kit LF00908U, LF00909U, LF00910U Instruction Manual

---





### Infant Circumcision Training Kit

- **LF00908U** White
- **LF00909U** Black
- **LF00910U** Brown

***All kits containing instruments are to be used for trainer purposes only. Instruments are not sterile and not intended for human use.***

## About the Simulator

With assistance from one of the top medical schools in Wisconsin, as well as medical professionals in South Africa and Indonesia, the **Life/form®** Infant Circumcision Trainer has been developed to offer realistic training of this delicate procedure. Now medical students and physicians can practice and improve their skills on the surgical removal of the foreskin on the glans penis without the worry of learning on a real live infant. This trainer may also be used to practice aftercare skills by the nursing staff and parents, and could be used to demonstrate the procedure to parents and family members. The skin is made from our lifelike material and is delicate and realistic to the touch. The trainer's foreskin is inexpensive and simple to replace.

### List of Components

- 6 glans with foreskins
- Infant trainer
- Nylon suture
- Dissection case
- Scalpel
- Forceps, micro
- Probe
- Scissors, Iris research
- Forceps, straight
- Forceps, 2 curved
- Blue soft case

### Available Supplies

- LF00904(A)U** White replacement glans & foreskin, pkg. of 10  
**LF00904(B)U** White replacement glans & foreskin, pkg. of 50  
**LF00900(A)U** Black replacement glans & foreskin, pkg. of 10  
**LF00900(B)U** Black replacement glans & foreskin, pkg. of 50  
**LF00902(A)U** Brown replacement glans & foreskin, pkg. of 10  
**LF00902(B)U** Brown replacement glans & foreskin, pkg. of 50  
**LF09919U** Nasco Cleaner

## General Instructions for Use

1. Place glans between thumb and forefinger and tug on foreskin to release skin from glans.
2. Press glans in hole in simulator.
3. Turn glans to be sure it is comfortably in place in simulator.
4. Tug slightly on glans to be sure it is tightly secure before procedure.
5. Remove glans after circumcision procedure.
6. Replace with another glans.

**Note:** Glans are not reusable. Replacement parts are sold in packages of 10 and 50.

## Care of the Simulator

Color pigments may vary per ethnic group. Normal surface soil can be removed with a mild soap and lukewarm water. Use Nasco Cleaner (LF09919U) to remove stubborn stains. Simply apply the cleaner to the soiled area and wipe with a clean, soft paper towel. **DO NOT** clean the simulator with solvents or corrosive materials, as they will cause damage.

## Other Available *Life/form*® Simulators

- LF00698U** Adult Injectable Arm (White)  
**LF00855U** Male Catheterization  
**LF00856U** Female Catheterization  
**LF00901U** Prostate Examination  
**LF00906U** Ostomy Care  
**LF00929U** Surgical Bandaging  
**LF00957U** Enema Administration  
**LF00958U** Pediatric Injectable Arm  
**LF00961U** Intramuscular Injection  
**LF00984U** Breast Examination  
**LF00995U** Arterial Puncture Arm  
**LF00999U** Pediatric Injectable Head  
**LF01005U** First Aid Arm  
**LF01008U** Intradermal Injection Arm  
**LF01012U** Heart Catheterization (TPN)  
**LF01019U** Ear Examination  
**LF01027U** Peritoneal Dialysis  
**LF01028U** Suture Practice Arm  
**LF01034U** Suture Practice Leg  
**LF01036U** Spinal Injection  
**LF01037U** Hemodialysis Practice Arm  
**LF01038U** Episiotomy Suturing Set  
**LF01042U** Suture Kit  
**LF01062U** Pelvic, Normal & Abnormal  
**LF01063U** Stump Bandaging, Upper  
**LF01064U** Stump Bandaging, Lower  
**LF01069U** Cervical Effacement  
**LF01070U** Birthing Station  
**LF01082U** Cricothyrotomy  
**LF01083U** Tracheostomy Care  
**LF01084U** Sigmoidoscopic Examination  
  
**LF01087U** Central Venous Cannulation  
**LF01095U** Blood Pressure Arm  
**LF01108U** Infant Intraosseous Infusion  
**LF01121U** Advanced IV Arm  
**LF01131U** Venipuncture and Injection Arm  
**LF01139U** Advanced IV Hand  
**LF01142U** Auscultation Trainer  
**LF01143U** Testicular Exam  
**LF01152U** Male & Female Catheter  
**LF01155U** Advanced CPR Dog  
**LF01162U** Venatech IV Trainer  
**LF01174U** NG Tube & Trach Skills  
  
**LF01184U** Venatech IM & Sub Q  
**LF01193U** Special Needs Baby  
**LF03000U** **CPARLENE**® Series  
**LF03601U** Adult Airway Management Trainer with Stand  
**LF03602U** Adult Airway Management Manikin  
**LF03609U** Child Airway Management Trainer with Stand  
**LF03616U** Child **CRiSis**™ Manikin  
**LF03617U** Deluxe Child **CRiSis**™ Manikin with Arrhythmia Tutor  
**LF03620U** PALS Update Kit  
**LF03623U** Infant Airway Management Trainer with stand  
**LF03632U** Child Intraosseous Infusion/ Femoral Access Leg on a Stand  
**LF03633U** Child Airway Management Trainer Torso  
**LF03693U** **Basic Buddy**® CPR Manikin  
**LF03699U** "Airway Larry" Airway Management Trainer  
**LF03709U** Infant **CRiSis**™ Manikin  
**LF03720U** **Baby Buddy**™ Infant CPR Manikin  
**LF03750U** Fat Old Fred  
**LF03760U** Airway Management/Cricoid Pressure Trainer  
  
**LF03770U** Chest Tube  
**LF03953U** **CRiSis**™ Manikin, complete  
**LF03955U** Deluxe **CRiSis**™ Manikin  
**LF03956U** Deluxe "Plus" **CRiSis**™ Manikin  
**LF03965U** Adult **CRiSis**™ Auscultation Manikin  
**LF03966U** Adult **CRiSis**™ Auscultation Manikin with ECG Simulator  
**LF04000U** **GERI**™/**KERI**™ Manikin Series  
**LF04200U** Adult Sternal Intraosseous Infusion  
**LF06001U** **CPR Prompt**® Adult/Child Manikin  
**LF06012U** **CPR Prompt**® Infant Manikin  
**LF06200U** **CPR Prompt**® Keychain Rescue Aid  
**LF06204U** **CPR Prompt**® Rescue and Practice Aid

*Nasco* Fort Atkinson

901 Janesville Avenue, P.O. Box 901  
Fort Atkinson, Wisconsin 53538-0901

1-800-558-9595

eNasco.com • E-mail: [lifeform@eNasco.com](mailto:lifeform@eNasco.com)

COPYRIGHT © NASCO 2010 PRINTED IN U.S.A. NP 062-11/RV 10-12